



## Shoe Size Charts

The Anatomical Concepts, Inc. line of orthoses has a number of adjustable features, two of which are adjustability for height and adjustability for length.

The chart below shows the adjustability limits for both height and length on each of the three sizes offered.

Shoe Size (USA)	Men's Heel to Toe (in inches)	Women's Heel to Toe (in inches)	Children Heel to Toe (in inches)
0			
1			3 9/16
2			3 7/8
3		8	4 1/4
3.5	8 3/4	8 3/16	4 7/16
4	8 15/16	8 5/16	4 9/16
4.5	9 1/16	8 1/2	4 11/16
5	9 1/4	8 11/16	4 7/8
5.5	9 5/16	8 13/16	5 1/16
6	9 3/8	9	5 1/4
6.5	9 1/2	9 3/16	5 3/8
7	9 11/16	9 5/16	5 9/16
7.5	9 13/16	9 1/2	5 3/4
8	10	9 11/16	5 7/8
8.5	10 3/16	9 13/16	6 1/16
9	10 5/16	10	6 1/4
9.5	10 1/2	10 3/16	6 3/8
10	10 11/16	10 5/16	6 9/16
10.5	10 13/16	10 1/2	6 3/4
11	11	10 11/16	6 7/8
11.5	11 3/16	10 13/16	7 1/16
12	11 5/16	11	7 1/4
12.5	11 1/2	11 3/16	7 3/8
13	11 11/16	11 5/16	7 9/16
13.5	11 13/16	11 1/2	7 3/4
14	12	11 11/16	
14.5	12 3/16		
15	12 5/16		
Size	Adult PRAFO	Pediatric PRAFO	Infant PRAFO
Height (in inches)	12 1/2 - 13 1/2	9 - 10	6 3/8 - 7 3/8
Length (in inches)	8 - 12 1/4	6 - 8 3/8	3 1/2 - 6 1/2