

The KMO™ (Knee Management Orthosis) features a single posterior joint that is a cost-effective, patient-friendly and a simple approach for static progressive positioning management of the knee.

Features

- Single posterior joint design allows for quick fit
- Patient friendly terry cloth liner system
- Pre and post operative management of bony and soft tissue injuries
- Adjustable range of stabilization of the knee in the sagittal and/or the coronal plane
- Static progressive positioning from 0° to 110°(additional ROM may be attained depending on size of patient's leg)
- Adjustable length from 12.50" to 16.50"(Pediatric) and from 17" to 26"(Adult)
- Available in Adult and Pediatric models, along with additional KMO liner variations



Other KMO™ Variations



**KMO™ (w/ PLASTIC THIGH & CALF SHELL)
KNEE ORTHOSIS**

**KMO™ EZ-FIT SYSTEM
KNEE ORTHOSIS**

Uses

- Chronic systemic conditions
- Immediate prophylactic protection at end range until treatment plans are finalized
- Contracture management
- CVA, TBI, SCI, CP, and MS
- Emergent lower limb immobilization
- Musculoskeletal injuries
- BKA(Below the Knee Amputation)
- Post Botox or serial cast management
- TKA(Total Knee Arthroplasty)
- Post-operative immobilization