

It's more than just positioning and offloading!

TERRA™

The TERRA™ Articulating Ankle/Foot Orthosis delivers similar benefits as the RIBBY™ Orthosis and ABBY™ Articulating Ankle/Foot Orthoses, but offers additional features appropriate for some patients. The integrity of the superstructure is unmatched by other Ankle/Foot orthoses



Features

- Clinically proven for effective ambulation*
- New non skid low profile plantar surface
- Open toe shoe provides at no extra cost
- Compatible with all liner variations that we offer if applicable
- Proven superior integrity of the superstructure**
- Available in Adult heavy duty #682NPB (Fits Women's shoe size range 5-14 and Mens shoe sizes 3.5-15)



Uses

- Rehab patients in orthotic care.
- Early intervention for the evolving rehab patient
 - CVA, TBA, SCI, MS, CP and Diabetic neuropathy
- Facilitation of functional knee and ankle joint position
 - Allows for optimal dorsiflexion ankle positioning, + or - 20° and Plantarflexion + or - 40° to provide adequate clearance during swing phase and knee stability in stance phase.
- Post musculoskeletal injury of the foot and ankle
 - Post TAL immobilization
 - Post Botox or serial casting
 - Stable fractures of the foot and ankle
 - Tendon or ligamentous tear/rupture

**DeToro W, (2001) Plantarflexion Resistance of Selected Ankle-Foot Orthoses: A Study of Commonly Prescribed Prefabricated and Custom-Molded Alternatives: J Prosthet Orthot:13;39-44

PDAC Approved L-1971

APU® Orthosis is a registered trademark of Anatomical Concepts Inc. #2,490,261

